

## Nursing Bench 2000



The nursing bench 2000 has a low starting height for easy transfer onto the bench as well as a perforated stretcher, which enables the water to pass through to the collection tray underneath.

The flexibility of the height adjustable nursing bench gives advantages to both the user and the carer, as the correct working height is equally important for the carer's safety and the user's comfort.




Electric with wired hand control – the red buttons clearly indicate "Up" and "Down". The wired hand control comes with a holder for fitting to the wall.



A water collection tray with integrated water outlet prevents water from splashing on carers or the floor.

The flexible hose attached to the centre mounted outlet of the water collection tray, can be connected to either an outlet in the floor or on the wall.

The integrated gas cylinders enable the bench to be folded out of the way when not in use to maximise floor space in front of the unit.



The stretcher of the bench is very easy to clean, and the stretcher cover can be unzipped for washing or easy replacement.

**Materials:** The frame is made from powder-lacquered stainless steel, the laying surface is a mesh cover. The water collection tray is made from polystyrene.

**Maximum load:** 175 kg/385 lbs

**Colour:** 000 White





### Nursing bench - adjustable height

- ◀ **R8468000-93** Electric nursing bench, 1800 mm/70"
- R8478000-93** Electric nursing bench with safety rail, 1800 mm/70"

When mounted at the recommended height, the height of the nursing bench can be adjusted from 300 mm/12" to 1000 mm/40" mm above floor.

If safety rail is to be pre-mounted, please be sure to indicate the correct item number when ordering.

Incl. wired hand control and flexible hose, 1500 mm/60".



### Safety rail










- RB2290** Length 1600 mm/63", for retro-fitting on nursing bench R8568000-93.

Provides safety for the user.

Can easily be folded up or down.



# Dimensioned drawings

Nursing Bench 3000		
<p>▼ R8520021 R8516021 R8539021 R8534021</p> 	<p>R8508021 R8503021 R8528021 ▼ R8523021</p> 	<p>R8510</p>  <p>R8511</p> 
Nursing Bench 2000		
<p>▼ R8468000-93 R8478000-93</p> 	<p>RB2290</p> 	
Nursing Bench 1000		
<p>▼ R8400000-93 R8402000-93</p> 	<p>R8401000-93 ▼ R8403000-93</p> 	<p>R8410000</p> 


Further information and mounting instructions can be found at [www.pressalit.com](http://www.pressalit.com)






# Mounting and accessories

## Mounting




R840000-93  
R8401000-93  
R8402000-93  
R8403000-93  
R8468000-93  
R8478000-93  
R8516021  
R8520021  
R8534021  
R8539021



R8503021  
R8508021  
R8523021  
R8528021


## Accessories

### RT610 / RT611



Shower head  
RT610: One spray mode  
RT611: Three spray mode

### RT640 / RT641



Shower mixer  
RT640: With thermostat  
RT641: With thermostat and extra long handles.

### RT018 / RT024 / RT036




Grab bar  
RT018xxx: 457,2 mm/18"  
RT024xxx: 610,0 mm/24"  
RT036xxx: 914,4 mm/36"  
Available in: white, anthracite, red and dark blue

### RT600



Showerhead holder  
For vertical handrail: RT042/RT118/  
RT119/RT136/RT138  
Available in: white, anthracite, red and dark blue

### R8420000 / R8421000



Splash curtain  
R8420000: Suitable for nursing bench R8400000-93 and R8402000-93  
R8421000: Suitable for nursing bench R8401000-93 and R8403000-93

Further information and mounting instructions can be found at [www.pressalit.com](http://www.pressalit.com)

